

Compensation Theory

Compensation control requirements for the proper instrument setup:

- **Unstained control** sample should be the **same origin** and prepared with the **same protocol** as the cells used in the experiment
- **Single-color controls:**
 - Need to have single-color controls for **all** the fluorochromes used in the experiment
 - All controls should be treated in the same way as samples (such as: incubation, fixation, permeabilization)
 - Should contain particles (either cells or compensation beads) stained with **each fluorescence marker individually**;
- **Note: Most common fluorescence markers:** *fluorochrome-labeled antibody (directly conjugated or via biotin-SA-fluorochrome system), fluorescent proteins, functional probes, etc.*
- Within each single-color control tube the **unstained and positive particles** have to be of the **same origin** (either cells or beads);
- **Level of fluorescence** in the controls should be **same or brighter** than in the samples.

Compensation goals — **unstained and stained** particles should show the **same mean value (or median for digital instruments) in all the "bystander" channels except the one "dedicated" for it.**

Bead Compensation Controls

Reason to use BD FACSComp Bead as Single Stain Controls:

- Cell numbers are **limiting**
- Antigen is **dimly** expressed
- Antibody has **low affinity** to receptor
- **Positive** population is **rare**
- Use **tandem** dyes

FACSComp Bead — Kits e.g. from BD, Biolegend, eBiosciences

Details:

- Each Kit contains **two vials** with bead polystyrene particles;
 - **Beads uncoated** that have no binding capacity; used for negative population
 - **Beads coated** with Goat **anti-Ig kappa**, which bind light chain-bearing immunoglobulin; used for positive population. Be sure to use same Fluorochrome-conjugated monoclonal antibodies as actually used for cell staining in your experiment
- **Note:** there are various FACSComp Bead Kits available with the coated beads of different specificity (Hamster, Rat or Mouse)

BD FACSComp Beads labeling procedure

(implemented and recommended by FCF-Berg):

- 1 **Verify origin** of your antibody (Hamster, Rat or Mouse) to specify which beads should be used
 - 2 **Mix Beads** well by shaking the bottle
 - 3 **Pre-dilute** anti-Ig Beads using 400 µl of PBS for 1 drop of Beads
 - 4 **Dispense** 100 µl of pre-diluted Beads per single color control sample
 - 5 **Add** 5 µl of pre-diluted **primary Ab** to the Beads
 - 6 **Incubate** 20 min at RT, mix occasionally
 - 7 If applicable, add 5 µl of pre-diluted secondary Ab with Fluorochrome
 - 8 **Incubate** 20 min at RT, mix occasionally
 - 9 **Add** 50 µl of pre-diluted **uncoated Beads** to each sample
 - 10 **Make** separate tube with **uncoated** beads for **Negative Control**
- Use** immediately for Multi-color compensation **or store** at +4°C for up to 2 weeks